


Research for lncRNA Library Generation kits

TaKaRa _ SMARTer pico

Hu Xi
2018-12-28

Mechanism of remove ribosomal cDNA(SMARTer pico kit)



- Removal of ribosomal cDNA after cDNA synthesis using probes specific to mammalian rRNA.
- These R-Probes hybridize to ribosomal RNA and mitochondrial rRNA sequences
- The ribosomal cDNA is then cleaved by **ZapR** in the presence of **the mammalian-specific R-Probes**

Contents

1. Summary of current methods for removing rRNA (principle)
 2. Summary and **comparison with SMARTer pico** kit
 3. Compare the **principle** of the two methods in the cDNA phase (PDD&DSN / DASH&Cas9)
 4. Probe **design** method for PDD&DSN
 5. Probe **design & generation** method for DASH&Cas9
 6. Compare the differences between the two methods
 7. Designed a feasible **test method**
- With Hongke Wang

1.Summary of current methods for removing rRNA (principle)

Depletion of abundant sequences

Summary from papers

Biotinylated capture probe — rRNA

RNase H & DNase I digestion — rRNA

Commercial kits


Illumina' s RiboZero/Qiagen GeneRead rRNA depletion
/Lexogen RiboCop /NEBNext rRNA depletion
/Roche /Kapa RiboErase /Takara/Clontech' s RiboGone

Probe-Directed Degradation (PDD)&Duplex-specific nuclease(DSN)

Depletion of Abundant Sequences by Hybridization (DASH)& CRISPR/Cas9


Biotinylated capture probe

Illumina RiboZero :




Thermo RiboMinus :

Step 1



Step 2




Step 3




RNase H & DNase I digestion

NEBNext rRNA Depletion Kit



Vazyme Ribo-off rRNA Depletion Kit

1. rRNA与探针结合 (rRNA probe hybridization)



2. RNase H消化 (RNase H depletion)




3. DNase I消化 (DNase I depletion)



- RNase H is an endoribonuclease that **specifically hydrolyzes the phosphodiester bonds of RNA which is hybridized to DNA**. This enzyme does not digest single or double-stranded DNA.



Probe-Directed Degradation (PDD)&Duplex-specific nuclease




Archer et al. BMC Genomics 2014

DSN(Duplex-specific nuclease)shows a strong preference for cleaving double-stranded (**ds**) **DNA** and **DNA in DNA-RNA hybrid duplexes**, and is practically inactive towards single-stranded (ss) DNA or single or double-stranded RNA.

Depletion of abundant sequences _ CRISPR/Cas9



Gu et al. *Genome Biology* 2016



Adapter-Dimer sgRNAs (targets library (-) strand)

20 bp Adapter Dimer sgRNAs
(AD1-AD4)

...TCCTACACGACGCTCTT **CCGATCTXXXXXACGGGCTAAT** ATTATCGGTGG...

5' Multiplex Adapter (X6 = 4 Unique 6 bp Sequences)

3' Adapter

hsa-miR-16-5p sgRNA (targets library (+) strand)

21 bp hsa-miR-16-5p
sgRNA

...TCCTACACGACGCTCTTCCGATCTXXXXXTA **GCAGCACGTAAATATTGGCGACGG** GCTAATATTATCGGTGG...

5' Multiplex Adapter

hsa-miR-16-5p

3' Adapter

NGG PAM


gRNA target sequence

Hardigan et al. *bioRxiv* 2018

2.Summary and comparison with SMARTer pico kit


Method	Stage	Input	Performance	Note
Biotinylated capture probe	RNA	>100ng	~ 87%	Biotinylated probe : rRNA / Streptavidin Beads Ribo-Zero™ Magnetic Gold Kit(Human/Mouse/Rat), 24 Rxns
RNase H & DNase I digestion	RNA	5ng to 1ug	~ 99%	Probes : rRNA / RNase H NEBNext rRNA Depletion Kit (Human/Mouse/Rat) 6 Rxns
PDD & DSN (2014)	ssDNA	500ng (start)	~ 94%	Probes : cDNA / DSN ; Incubation temperature : 75°C—48°C (Gradient cooling) , 1h Termination condition : stop solution
DASH&Cas9 *	dsDNA	<10ng	~ 99%	sgRNA : cDNA / “NGG” PAM motif : Cas9 Incubation temperature : 37°C , 2h ; Termination condition : heat or stop solution
TaKaRa/Clontech SMARTer pico**	dsDNA	250 pg to 10 ng	~ 95%	R-Probes : cDNA / ZapR Incubation temperature : 37°C , 1h ; Termination condition : heat

*samples from plasma ** E-mail consultation





Roche NimbleGen SeqCap EZ Library SR

Both the Exome and Clontech methods were effective and sensitivity.



Zachary T. Herbert et al. BMC Genomics 2018

3. Compare the **principle** of the two methods in the cDNA phase (PDD&DSN / DASH&Cas9)



INPUT

cDNA Library form
Reverse Transcription



PDD&DSN

RNase A remove
RNA from cDNA
library



Hybrid incubation
and DSN treatment



Purification of cDNA



DASH&Cas9

Hybrid incubation
with sgRNA and Cas9



RNase A remove
sgRNA



Purification of cDNA

OUTPUT



PCR Amplification



Gel separation
selection wanted size
sequence

4.PDD&DSN

Probes design

NCBI BLAST/ ApE (A Plasmid Editor, v2.0.47)/ Biopython (www.biopython.org)

Define number of PDD probes

Calculate optimal probe density

The maximum inter-probe distance = the minimum library insert size - 20 nt

The period (distance from the start of one probe to the start of the next) = this distance + the minimum probe length

Allow coverage of nearly 10 kb of target sequence using 50 probes

Identify regions of cross-homology in the target transcript



(www.ncbi.nlm.nih.gov/blast)

Enter
the target sequence

Database: refseq_rna

Choose the species

“Optimize for” → “somewhat similar sequences (blastn)”

Algorithm parameters:

“Max target sequences” to 10000;
“Expect threshold” to 100000;
“word size” to 7;
“match/mismatch scores” to 1/-4;
“Gap costs” to “existence: 5 extension: 2”

BLAST

Results export
as an “XML”

find short stretches of high homology to other transcripts in the target RNA

Highlight cross-homology in ApE

Python script: [blast2ape.py](https://github.com/stu2/xml2ape) <https://github.com/stu2/xml2ape>

In ApE:

Paste the target RNA sequence → save it

Download the BLAST hit table text file

Add the xml2ape.Py script

Execute the script

Produce an ApE file

Select probes in ApE


5.DASH&Cas9

sgRNAs design & generation

<https://zlab.bio/guide-design-resources>



SnapGene
display design rules



Find and identify targeting sequences



Design sgRNA Transcription template



Synthetic DNA template



PCR template amplification




In vitro transcription (IVT)



Digesting DNA



Purified RNA



6. Compare the differences between the two methods

CONTENT	PDD&DSN	DASH & Cas9
Probes types	DNA, ~15nt	RNA, ~20nt (targets sequences)
Probes preparation method	Company synthesis	Company synthesis & Self preparation
Probe storage condition	-20C°, 12 months	-80C°, 6 months
Target sequence sequence condition	>35nt, Suitable for lncRNA-seq	no request, Suitable for lncRNA-seq & smRNA-seq
Ability to remove Adapter dimer**	YES	YES
Key reagents	DSN(Evrogen, EA003-10U): ¥1820/10 reactions HEPES-KOH (Sigma, H0527-25G): ¥644.67	MEGAscript™ T7 Transcription Kit (Invitrogen, AM1333): ¥3122/25 reactions TURBO™ DNase (Invitrogen, AM2238): ¥1805/500 reactions MEGAclean™ Transcription Clean-Up Kit (Invitrogen, AM1908): ¥1720/50 reactions S. pyogenes Cas9(PNABio, CPO2-250ug): ¥2540.27/50 reactions NEBuffer™ 3.1(B7203S): ¥235.2/2500 reactions

**Adapter dimer vs Primer dimer

~ ¥10/probe

~ ¥182/reaction

~ ¥220 (~ ¥85_TaKaRa) /probe*

~ ¥52/reaction

*idea : Design universal primers to amplify DNA templates

7. Designed a feasible test method

Samples : Standard RNA, such as **Spike In**

Methods : Choose two RNA sequences

rS1=Control sequences

rS2=Target sequences

Probes design

Reverse Transcription / RT-PCR

DSN / Cas9

PCR Amplification/qPCR

Gel electrophoresis/Quantitative analysis

